

Quant-iT™ Assays

Abbreviated Protocol

NOTE: For best results, store the dye and the buffer at room temperature. Store the DNA, RNA and protein standards at 4°C. Ensure that all assay reagents are at **room temperature** before you begin.

1. Set up two tubes for the standards (three for the protein assay) and one tube for each user sample.
2. Prepare the Quant-iT™ **Working Solution** by diluting the Quant-iT™ reagent 1:200 in Quant-iT™ buffer. 200 µL of **Working Solution** is required for each sample and standard.
3. Prepare the Assay Tubes* according to the table below.


	Standard Assay Tubes	User Sample Assay Tubes
Volume of Working Solution (from step 2) to add	190 µL	180–199 µL
Volume of Standard (from kit) to add	10 µL	—
Volume of User Sample to add	—	1–20 µL
Total Volume in each Assay Tube	200 µL	200 µL

4. Vortex all tubes for 2–3 seconds.
5. Incubate the tubes for 2 minutes at room temperature (15 minutes for the Quant-iT™ protein assay).
6. Insert the tubes in the Qubit® fluorometer and take readings.
7. Multiply the reading by the dilution factor (see Manual) to determine concentration of your original sample. Alternatively, choose **Calculate sample concentration** to have the Qubit® fluorometer perform this multiplication for you.

* Use only thin-wall, clear 0.5 mL PCR tubes. Acceptable tubes include Qubit™ assay tubes (500, Invitrogen Cat. no. Q32856) or Axygen PCR-05-C tubes (VWR, part number 10011-830).

www.invitrogen.com/qubit

Ensure all reagents are at room temperature



* where n = number of Standards plus number of Samples

 invitrogen™

Vortex all assay tubes for 2–3 seconds

Incubate at room temperature for 2 minutes (15 minutes for Quant-iT™ protein assay)

Read tubes in Qubit® fluorometer

